



**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
East Avenue, Quezon City 1100  
981-0300 / 981-0400 local 1157  
<http://www.nkti.gov.ph/>

**SUPPLEMENTAL BID BULLETIN NO. 22-013-1**

**SUPPLY AND DELIVERY OF VARIOUS RADIOLOGY SUPPLIES (LINE BIDDING)**

**Bid Reference IB No. 22-013**


This Supplemental Bid Bulletin No. 22-013-1 is being issued to clarify, modify and amend items/specifications in the Bid Documents for the aforementioned project.

- a. **Section VII Revised Technical Specifications**  
*(Please see attached file)*
- b. **Revised Price Schedule "Annex A"**  
*(Please see attached file)*

This Supplemental Bid Bulletin including Annexes, if any, shall form part of the Bid Documents. Any provisions in the Bid Documents inconsistent herewith is hereby amended, modified and superseded accordingly.

For guidance and information of all concerned.

Issued this 28<sup>th</sup> day of October 2021 in Quezon City.

  
**ARNOLD JOSEPH M. FERNANDEZ, MD**  
**BAC, Chairman**

Received by: (PLS SIGN) \_\_\_\_\_

Bidder's Name: (PLS PRINT) \_\_\_\_\_

Date: \_\_\_\_\_



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
East Avenue, Quezon City 1100  
8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

## ***Section VII.***

### ***Revised Technical Specifications***

***NKTI Reference No. IB 22-013***



## Technical Specifications

**Instruction:**

*Bidders must state in the column provided either “Comply” or “Not Comply” against each of the individual parameters of each Specification stating the corresponding performance parameter of the equipment offered. Statements of “Comply” or “Not Comply” must be supported by evidence in a Bidders Bid and cross-referenced to that evidence. Evidence shall be in the form of manufacturer’s un-amended sales literature, unconditional statements of specification and compliance issued by the manufacturer, samples, independent test data etc., as appropriate. A statement that is not supported by evidence or is subsequently found to be contradicted by the evidence presented will render the Bid under evaluation liable for rejection. A statement either in the Bidder’s statement of compliance or the supporting evidence that is found to be false either during Bid evaluation, post-qualification or the execution of the Contract may be regarded as fraudulent and render the Bidder or supplier liable for prosecution subject to the applicable laws and issuances.*

Item No.	Particulars	Specification	Statement of Compliance
1	Micro-Catheter, Angiographic Fr. 2.4 – 3.0 with Guidewire (Co-axial Catheter System)	<p><b>Packaging:</b></p> <p><b>Content:</b></p> <p>a) Microcatheter</p> <p>b) Guidewire</p> <p>-Single use, disposable</p> <p>-Sterility: Sterile</p> <p>-Individual, see through with peel away pack</p> <p>-Label: Lot/batch number, expiration date</p> <p><b>Specifications:</b></p> <p><b>a)Microcatheter:</b></p> <p>-130 cm in length</p> <p>-PTFE Platinum/ iridium marker</p> <p>-Hydrophilic coating on the entire surface over its entire length except its proximal end coating that gives lubricity when its wet</p> <p><b>b) Guidewire:</b></p> <p>-140-165 cm in length</p> <p>-Nitinol, hydrophilic coating on its surface and adequate length of radiopaque tip</p> <p>-Preshaped type/ Shapeable tip</p> <p><b>Functionality:</b></p> <p><b>-Durability:</b> Maintain stability through out the procedure</p> <p style="text-align: right;">: Not easily kinked or broken</p>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		<ul style="list-style-type: none"> <li>-Should be compatible with lipiodol, alcohol glue, dimethyl sulfoxide, various coils, embolic spheres, radiomimetic substances and radiation spheres</li> <li>-Guidewire should have frictionless delivery of coils and other embolics, distal length should be shapeable.</li> </ul>	
2	Catheter, All purpose drainage, Single F. 12	<p><b>Packaging:</b></p> <p><b>CONTENT: (catheter, puncture needle and connecting tube)</b></p> <ul style="list-style-type: none"> <li>-Sterile</li> <li>-Individual, see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-At least 25-30 cm working length</li> <li>-Radiopaque shaft, hydrophilic coating</li> <li>-Accommodate guidewire size 0.035-0.038 inches</li> <li>-Distal end contains a self-retaining “pigtail”</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-No links, breakage and maintain stability throughout the procedure</li> <li>-Catheter tip should be atraumatic, smooth, firm but flexible</li> <li>-Should have a direct access (single stick method) type</li> <li>-Should have a good locking mechanism to fasten the catheter and prevent accidental pull-out</li> <li>-Distal end contains a self-retaining “pigtail”</li> </ul>	
3	Catheter, All purpose drainage, Set F. 12	<p><b>KIT/SET CONTENT:</b></p> <ul style="list-style-type: none"> <li>-1 piece trocar needle</li> <li>-1 piece stiff guidewire 0.038 inch</li> <li>-1 piece multi-purpose drainage catheter</li> <li>-1 piece universal connecting tube (from catheter to any collecting bag)</li> <li>-1 piece stopcock</li> <li>-3 dilators beginning with Fr. 6, 8 and 10</li> </ul> <p><b>Packaging:</b></p>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		<ul style="list-style-type: none"> <li>-Sterility: Sterile</li> <li>-See through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-25-30 cm working length</li> <li>-Radiopaque shaft, hydrophilic coating catheter</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-No kinks, breakage and maintain stability throughout the procedure</li> <li>-Catheter tip should be atraumatic, smooth, firm but flexible</li> <li>-Should have a direct access (single stick method) type</li> </ul>	
4	Catheter, Angiographic Cobra 2 Fr. 5 (middle curve) (C2)	<p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Single use, disposable</li> <li>-Sterility: Sterile</li> <li>-Individual, see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-At least 80 cm in length</li> <li>-No side holes</li> <li>-Open end hole</li> <li>-Polymer, stainless steel braided wire</li> <li>-Hydrophilic, radiopaque and torqueable catheter</li> <li>-Luer lock hub at the proximal end</li> <li>-Accommodate up to 0.038 inch guidewire size</li> <li>-Up to 1000 PSI</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-<b>Durability:</b> Maintain stability throughout the procedure : Not easily kinked or bent</li> <li>-<b>Ease of use:</b> Atraumatic and smooth insertion</li> </ul>	
5	Catheter, Angiographic Yashiro F. 5	<p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Single use, disposable</li> <li>-Sterile</li> <li>-Individual, see through with peel away pack</li> </ul>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		<p>-Label: Lot/batch number, expiration date</p> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-at least 65-80 cm in length</li> <li>-No side holes</li> <li>-Open end hole</li> <li>-Polymer, Stainless steel braided wire</li> <li>-Hydrophilic, radiopaque and torqueable catheter</li> <li>-Luer lock hub at the proximal end</li> <li>-Accommodate up to 0.038 inch guidewire size</li> <li>-Up to 1000 PSI</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-<b>Durability:</b> Maintain stability throughout the procedure : Not easily kinked or bent</li> <li>-<b>Ease of use:</b> Atraumatic and smooth insertion</li> </ul>	
6	Guidewire, J-tip	<p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Usage: Single use (disposable), medical</li> <li>-Sterility: Sterile</li> <li>-Individual, double see through with peel away pack</li> <li>-Label: Batch/lot number, expiry date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-0.038 inch diameter</li> <li>-150 cm in length</li> <li>-Stainless steel, movable core</li> <li>-Hydrophilic coating / Nitinol</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-<b>Durability:</b> Kink resistant, peel resistant/ flexible guidewire</li> <li>-<b>Ease of use:</b> Allows smooth passage and exchange between catheters, cannula and tissue</li> <li>-<b>Identification:</b> Radiopaque</li> </ul>	
7	Gelatin, Absorbable Hemostat Sponge	<p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Single use, disposable</li> <li>-Sterility: Sterile</li> <li>-Individual, double see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		<p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-80mm x 50mm x 10mm</li> <li>-Contain highly purified neutral gelatin foam of uniform fine porosity</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-<b>Ease of use:</b> Can be cut to any shape or sizes</li> <li>-Quickly adheres to the wound surface and control bleeding</li> </ul>	
8	Needle, IV Cannula G. 20	<p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Single use, disposable</li> <li>-Sterile</li> <li>-Individual, see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-32mm length</li> </ul> <p><b>Features:</b></p> <ul style="list-style-type: none"> <li>-With luer lock plug/cap</li> <li>-With safety feature against needleprick</li> <li>-With injection port for medicine and contrast media</li> <li>-Universal color-coding of needle-bore gauges</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-Flexible and durable catheter</li> <li>-<b>Ease of use:</b> Sharp needle, smooth and glides easily</li> </ul>	
9	Needle, IV Cannula G. 22	<p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Single use, disposable</li> <li>-Sterile</li> <li>-Individual, see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-25mm length</li> </ul> <p><b>Features:</b></p> <ul style="list-style-type: none"> <li>-With luer lock plug/cap</li> <li>-With safety feature against needleprick</li> <li>-With injection port for medicine and contrast media</li> <li>-Universal color-coding of needle-bore</li> </ul>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		gauges <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
10	Needle, IV Cannula G. 24	<b>Packaging:</b> -Single use, disposable -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -19mm length <b>Features:</b> -With luer lock plug/cap -With safety feature against needleprick -With injection port for medicine and contrast media -Universal color-coding of needle-bore gauges <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
11	Needle, Biopsy Chiba G. 22 X 20cm	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -20 cm in length (8") -Quincke-type bevel tip <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
12	Needle, Biopsy Chiba G. 22 x 10cm	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		-10 cm in length (4") -Quincke-type bevel <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
13	Needle, Biopsy Chiba G. 18 x 20cm	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -20 cm in length -Quincke-type bevel with echogenic marker -Stainless steel <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
14	Needle, Centesis Catheter F. 5	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -At least 10 cm in length <b>Features:</b> -Tapered catheter tip: Decrease resistance for less patient discomfort -Radiopaque / Echogenic needle tip: Visibility -Has a four radial side holes on the catheter – for maximum drainage -Clear needle hub: Visual flashback confirmation -Vent plug: Prevent fluid leakage <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
15	Needle Core Biopsy G. 16 (semi-automatic)	<b>Packaging:</b> -Sterile	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		<p>-Individual, see through with peel away pack</p> <p>-Label: Lot/batch number, expiration date</p> <p><b>Specifications:</b></p> <p>-15 cm in length</p> <p><b>Features:</b></p> <p>-Variable notch sizes, 10mm or 20mm allows for precise sampling for various bx procedures</p> <p>-Has a short four-sided needle bevel for quick easy insertion with less trauma</p> <p>-With echogenic marks on the stylet and cannula provide visibility under ultrasound</p> <p><b>Functionality:</b></p> <p><b>-Durability:</b> Flexible and durable catheter</p> <p><b>-Ease of use:</b> Sharp needle, smooth and glides easily</p>	
16	Needle, Coaxial Introducer G. 15	<p><b>Packaging:</b></p> <p>-Sterile</p> <p>-Individual, see through with peel away pack</p> <p>-Label: Lot/batch number, expiration date</p> <p><b>Specifications:</b></p> <p>-10 cm in length</p> <p>-Whitacre bevel tip/ pencil point</p> <p><b>Functionality:</b></p> <p>-Flexible and durable catheter</p> <p><b>-Ease of use:</b> Sharp needle, smooth and glides easily</p> <p>-Should be compatible in the core needle bx</p>	
17	Needle, Core Biopsy G. 18 (semi-automatic)	<p><b>Packaging:</b></p> <p>-Sterile</p> <p>-Individual, see through with peel away pack</p> <p>-Label: Lot/batch number, expiration date</p> <p><b>Specifications:</b></p> <p>-15 cm in length</p> <p><b>Features:</b></p> <p>-Variable notch sizes, 10mm or 20mm allows for precise sampling for various bx procedures</p> <p>-Has a short four-sided needle bevel for</p>	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		quick easy insertion with less trauma -With echogenic marks on the stylet and cannula provide visibility under ultrasound <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
18	Needle, Coaxial Introducer G. 17	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -10 cm in length -Whitacre bevel tip <b>Functionality:</b> -flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
19	Needle, Core Biopsy G. 20 (semi-automatic)	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -15 cm in length <b>Features:</b> -Variable notch sizes, 10mm or 20mm allows for precise sampling for various bx procedures -Has a short four-sided needle bevel for quick easy insertion with less trauma -With echogenic marks on the stylet and cannula provide visibility under ultrasound <b>Functionality:</b> -Flexible and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
20	Needle, Coaxial Introducer G. 19	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack	



**Republic of the Philippines**  
**National Kidney and Transplant Institute**  
**Bids and Awards Committee**  
 East Avenue, Quezon City 1100  
 8981-0300 / 9881-0400 local 1156/1186  
<http://www.nkti.gov.ph/>

		-Label: Lot/batch number, expiration date <b>Specifications:</b> -10 cm in length -Whitacre bevel tip <b>Functionality:</b> -Flexibe and durable catheter <b>-Ease of use:</b> Sharp needle, smooth and glides easily	
21	Disinfectant, wipes	<b>Packaging:</b> -Preferably in plastic canister -At least 160 sheets -Resealable or flip top cap (easy to open and sealed cap) <b>Specifications:</b> -Non-woven towelettes -Active ingredients: Di isobutyl phenoxy ethoxy ethyl dimethyl benzyl ammonium chloride <b>Functionality:</b> -Disinfection and decontamination of medical devices and radiology equipment/machine -Easy to store	

**We guarantee that we will submit sample within 3 calendar days upon notice from the BAC. Failure to submit within the prescribed period will automatically disqualify the bidder.**

*Lowest Calculated Bidder shall also submit clear copies of the colored picture and **COMPLETE PRODUCT DATA SHEET** [i.e brochures, pamphlets, sales literature, documentation of standards of compliance, etc.] with the required specifications already highlighted of the item proposed in A4 size photo paper indicating the brand name with reference no.*

***(NOTE: Submission of samples are no longer necessary if used by the end-user or otherwise specified by the end-user).***

<p><b>Conforme:</b></p> <p>Name: _____</p> <p>Legal Capacity: _____</p> <p>Signature: _____</p> <p>Duly authorized to sign the Bid for and behalf of: _____</p>
---

REVISED PRICE SCHEDULE

"ANNEX A"

IB No. 22-013: SUPPLY AND DELIVERY OF VARIOUS RADIOLOGY SUPPLIES (Line Bidding)

BAC Conference Room, G/F NKT Main Building, East Avenue, Quezon City

ITEM NO.	ITEM / SUPPLIES	TECHNICAL SPECIFICATION	UOM	TOTAL QTY.	NKT APPROVED BUDGET COST		BIDDER'S PROPOSAL		BRAND	REMARKS
					UNIT COST	TOTAL COST	UNIT COST	TOTAL COST		
1	Micro-Catheter, Angiographic Fr. 2.4-3.0 with Guidewire (Co-axial Catheter System)	<p><b>Packaging:</b></p> <p><b>Content:</b></p> <ul style="list-style-type: none"> <li>a) Microcatheter</li> <li>b) Guidewire</li> <li>-Single use, disposable</li> <li>-Sterility: Sterile</li> <li>-Individual, see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <p><b>a) Microcatheter:</b></p> <ul style="list-style-type: none"> <li>-130 cm in length</li> <li>-PTFE Platinum/ Iridium marker</li> <li>-Hydrophilic coating on the entire surface over its entire length except its proximal end coating that gives lubricity when its wet</li> </ul> <p><b>b) Guidewire:</b></p> <ul style="list-style-type: none"> <li>-140-165 cm in length</li> <li>-Nitinol, hydrophilic coating on its surface and adequate length of radiopaque tip</li> <li>-Preshaped type/ Shapeable tip</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-<b>Durability:</b> Maintain stability through out the procedure : Not easily kinked or broken</li> <li>-Should be compatible with lipiodol, alcohol glue, dimethyl sulfoxide, various coils, embolic spheres, radiomimetic substances and radiation spheres</li> <li>-Guidewire should have frictionless delivery of coils and other embolics, distal length should be shapeable.</li> </ul>	kit/set	120	14,300.00	1,716,000.00				
2	Catheter, All purpose drainage, Single F. 12	<p><b>Packaging:</b></p> <p><b>CONTENT: (catheter, puncture needle and connecting tube)</b></p> <ul style="list-style-type: none"> <li>-Sterile</li> <li>-Individual, see through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-At least 25-30 cm working length</li> <li>-Radiopaque shaft, hydrophilic coating</li> <li>-Accommodate guidewire size 0.035-0.038 inches</li> <li>-Distal end contains a self-retaining "pigtail"</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-No links, breakage and maintain stability throughout the procedure</li> <li>-Catheter tip should be atraumatic, smooth, firm but flexible</li> <li>-Should have a direct access (single stick method) type.</li> <li>-Should have a good locking mechanism to fasten the catheter and prevent accidental pull-out</li> <li>-Distal end contains a self-retaining "pigtail"</li> </ul>	pack	12	2,800.00	33,600.00				
3	Catheter, All purpose drainage, Set F. 12	<p><b>KIT/SET CONTENT:</b></p> <ul style="list-style-type: none"> <li>-1 piece trocar needle</li> <li>-1 piece stiff guidewire 0.038 inch</li> <li>-1 piece multi-purpose drainage catheter</li> <li>-1 piece universal connecting tube (from catheter to any collecting bag)</li> <li>-1 piece stopcock</li> <li>-3 dilators beginning with Fr. 6, 8 and 10</li> </ul> <p><b>Packaging:</b></p> <ul style="list-style-type: none"> <li>-Sterility: Sterile</li> <li>-See through with peel away pack</li> <li>-Label: Lot/batch number, expiration date</li> </ul> <p><b>Specifications:</b></p> <ul style="list-style-type: none"> <li>-25-30 cm working length</li> <li>-Radiopaque shaft, hydrophilic coating catheter</li> </ul> <p><b>Functionality:</b></p> <ul style="list-style-type: none"> <li>-No kinks, breakage and maintain stability throughout the procedure</li> <li>-Catheter tip should be atraumatic, smooth, firm but flexible</li> <li>-Should have a direct access (single stick method) type</li> </ul>	set/kit	12	4,900.00	58,800.00				

4	Catheter, Angiographic Cobra 2 French 5 (middle Curve) (C2)	<b>Packaging:</b> -Single use, disposable -Sterility: Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -At least 80 cm in length -No side holes -Open end hole -Polymer, stainless steel braided wire -hydrophilic, radiopaque and torqueable catheter -Luer lock hub at the proximal end -Accommodate up to 0.038 inch guidewire size -Up to 1000 PSI <b>Functionality:</b> <b>-Durability:</b> Maintain stability throughout the procedure : Not easily kinked or bent <b>-Ease of use:</b> Atraumatic and smooth insertion	piece	120	2,500.00	300,000.00				
5	Catheter, Angiographic Yashiro F.5	<b>Packaging:</b> -Single use, disposable -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> <b>-at least 65-80</b> cm in length -No side holes -Open end hole -Polymer, Stainless steel braided wire -hydrophilic, radiopaque and torqueable catheter -Luer lock hub at the proximal end -Accommodate up to 0.038 inch guidewire size -Up to 1000 PSI <b>Functionality:</b> <b>-Durability:</b> Maintain stability throughout the procedure : Not easily kinked or bent <b>-Ease of use:</b> Atraumatic and smooth insertion	piece	120	2,500.00	300,000.00				
6	Guidewire, J-Tip	<b>PACKAGING:</b> -Usage: single use (disposable), medical -Sterility: sterile -Individual, double see through with peel away pack -Label: batch/lot number, expiry date <b>SPECIFICATIONS:</b> -0.038 inch diameter -150cm in length -Stainless steel, movable core -hydrophilic coating / Nitinol <b>FUNCTIONALITY:</b> Durability - Kink resistant, peel resistant/ flexible guidewire Ease of Use - Allows smooth passage and exchange between catheters, cannula, and tissue Identification - Radiopaque	piece	120	1,000.00	120,000.00				
7	Gelatin, Absorbable Hemostat Sponge	<b>Packaging:</b> -Single use, disposable -Sterility: Sterile -Individual, double see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -80mm x 50mm x 10mm -Contain highly purified neutral gelatin foam of uniform fine porosity <b>Functionality:</b> -Ease of use: Can be cut to any shape or sizes -Quickly adheres to the wound surface and control bleeding	pack	180	380.00	68,400.00				

8	NEEDLE, IV CANNULA G.20	<b>PACKAGING:</b> -Single use, disposable. -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>SPECIFICATIONS:</b> -32mm length <b>FEATURES:</b> -With luer lock plug/cap -With safety feature against needleprick -With injection port for medicine and contrast media -Universal color-coding of needle-bore gauges <b>FUNCTIONALITY:</b> -Flexible and durable catheter -Ease of use: sharp needle, smooth and glides easily	piece	4000	80.00	320,000.00			
9	Needle, IV Cannula G. 22	<b>Packaging:</b> -Single use, disposable -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -25mm length <b>Features:</b> -With luer lock plug/cap -With safety feature against needleprick -With injection port for medicine and contrast media -Universal color-coding of needle-bore gauges <b>Functionality:</b> -Flexible and durable catheter -Ease of use: Sharp needle, smooth and glides easily	piece	800	80.00	64,000.00			
10	Needle, IV Cannula G. 24	<b>Packaging:</b> -Single use, disposable -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -19mm length <b>Features:</b> -With luer lock plug/cap -With safety feature against needleprick -With injection port for medicine and contrast media -Universal color-coding of needle-bore gauges <b>Functionality:</b> -Flexible and durable catheter -Ease of use: Sharp needle, smooth and glides easily	piece	800	80.00	64,000.00			
11	NEEDLE, BIOPSY CHIBA G.22 x 20cm	<b>PACKAGING:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>SPECIFICATIONS:</b> -20cm in length (8") -Quincke-type bevel tip <b>FUNCTIONALITY:</b> -Flexible and durable catheter -Ease of use: sharp needle, smooth and glides easily	piece	240	750.00	180,000.00			
12	Needle, Biopsy Chiba G. 22 x 10cm	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -10 cm in length (4") -Quincke-type bevel <b>Functionality:</b> -Flexible and durable catheter -Ease of use: Sharp needle, smooth and glides easily	piece	36	1,000.00	36,000.00			

13	Needle, Biopsy Chiba G. 18 x 20 cm	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -20cm in length -Quincke-type bevel with echogenic marker -Stainless steel <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	240	1,000.00	240,000.00				
14	Needle, Centesis Catheter F. 5	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> - at least 10 cm in length <b>Features:</b> -Tapered catheter tip: decrease resistance for less patient discomfort. -Radiopaque / Echogenic needle tip: visibility -Has a four radial side holes on the catheter - for maximum drainage -Clear needle hub: visual flashback confirmation -Vent plug: Prevent fluid leakage <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	120	1,920.00	230,400.00				
15	Needle Core Biopsy G. 16 (semi-automatic)	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -15 cm in length <b>Features:</b> -Variable notch sizes, 10mm or 20mm allows for precise sampling for various bx procedures -Has a short four-sided needle bevel for quick easy insertion with less trauma -With echogenic marks on the stylet and cannula provide visibility under ultrasound <b>Functionality:</b> -Durability: Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	60	2,500.00	150,000.00				
16	Needle, Coaxial Introducer G. 15	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -10cm in length -Whitacre bevel tip/ pencil point <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily -Should be compatible in the core needle bx	piece	60	1,250.00	75,000.00				
17	Needle, Core Biopsy G. 18 (semi-automatic)	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -15 cm in length <b>Features:</b> -Variable notch sizes, 10mm or 20mm allows for precise sampling for various bx procedures -Has a short four-sided needle bevel for quick easy insertion with less trauma -With echogenic marks on the stylet and cannula provide visibility under ultrasound <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	150	2,500.00	375,000.00				

18	Needle, Coaxial Introducer G. 17	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -10cm in length -Whitacre bevel tip <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	150	1,250.00	187,500.00					
19	Needle, Core Biopsy G. 20 (semi-automatic)	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -15 cm in length <b>Features:</b> -Variable notch sizes, 10mm or 20mm allows for precise sampling for various bx procedures -Has a short four-sided needle bevel for quick easy insertion with less trauma -With echogenic marks on the stylet and cannula provide visibility under ultrasound <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	72	2,500.00	180,000.00					
20	Needle, Coaxial Introducer G. 19	<b>Packaging:</b> -Sterile -Individual, see through with peel away pack -Label: Lot/batch number, expiration date <b>Specifications:</b> -10 cm in length -Whitacre bevel tip <b>Functionality:</b> -Flexible and durable catheter <b>Ease of use:</b> Sharp needle, smooth and glides easily	piece	72	1,250.00	90,000.00					
21	Disinfectant, wipes	<b>Packaging:</b> -Preferably in plastic canister -At least 160 sheets -Resealable or flip top cap (easy to open and sealed cap) <b>Specifications:</b> -Non-woven towelettes -Active ingredients: Di isobutyl phenoxy ethoxy ethyl dimethy benzyl ammonium chloride <b>Functionality:</b> -Disinfection and decontamination of medical devices and radiology equipment/machine -Easy to store	CANISTER	200	960.00	192,000.00					
						<b>GRAND TOTAL</b>	<b>4,980,700.00</b>				

\_\_\_\_\_

Printed Name of Company

\_\_\_\_\_

Date

\_\_\_\_\_

Signature

\_\_\_\_\_

Address

\_\_\_\_\_

Tel. No.

\_\_\_\_\_

Printed Name and Designation